



## **LOW PROFILE**

- Optimized folding process for low balloon profile
- Exceptional entry profile for successful lesion crossing
- 6.0 x 200 mm & 7.0 x 80 mm through 5Fr, 10.0 x 80 mm through 6Fr

## **EXCELLENT DELIVERABILITY**

- Slender tapered & atraumatic tip design for higher flexibility
- Hydrophilic coating for enhanced trackability

## ULTIMATE SOLUTION FOR CHALLENGING LESIONS

- High-pressure balloon up to 20 atm
- Extended size matrix: balloon size up to 200 mm
- Shafts 40, 75 & 135 cm

## Efficiency at its best

